

# Press Pelletizing Starter Package

Creating an optimal pressed pellet for XRF analysis that meets with standard ASTM and ISO methods requires a laboratory grade press as well as premium quality accessories and consumable items.

## Manual Press



## Die Set



## Filling Tool




## LIMITED TIME

### Starter Package Includes:

- 25 Ton Manual Press
- Die Set
- Filling Tool
- Aluminum Cups
- Pellet Film
- Bead Labels
- PlumaSol
- Pellet Blend Binder

## PlumaSol

Specialty grinding solution additive to improve cleaning during particle size reduction.



## Pellet Film



## Pellet Blend

- Excellent choice for Iron/Ore, Slag and Magnesites
- Exhibits excellent binding characteristics
- Suitable for mixing or grinding inhomogeneous samples
- Suitable to be used as a binder, referenced sample or calibration standard.



## Aluminum Cups

- Finest quality
- Straight & Taper Wall Design
- Thin wall material ensures minimal back pressure
- Resistance to accidental damage
- Easy to handle, label, and store



Request a quote: [sales@premierlabsupply.com](mailto:sales@premierlabsupply.com)

691 NW Enterprise Drive, Port Saint Lucie, Florida, 34986, USA

T. 772-873-1700 F. 772-873-1800, [www.premierlabsupply.com](http://www.premierlabsupply.com)